

| Acids | | | | | | | | | | | |
|-------------------------|---|----------------|-------|----|----|----|----|----|-----|----|----|
| Acetic acid | CH ₃ COOH | Conc. > 98% | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 | |
| Chromic acid | H ₂ CrO ₄ | 40% | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A6 | B1 | C0 | D0 | E0 | F0 | G0 | |
| Citric acid | C ₆ H ₈ O ₇ | 50% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Hydrochloric acid | HCl | Conc. 37% | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A5 | B1 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A5 | B1 | C0 | D0 | E0 | F0 | G0 | |
| Hydrofluoric acid | HF | 40% | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A1 | B1 | C1 | D1 | E0 | F0 | G0 | |
| Phosphoric acid | H ₃ PO ₄ | Conc. > 85% | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B2 | C1 | D0 | E0 | F0 | G0 | |
| Phosphoric acid | H ₃ PO ₄ | 38% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Lactic acid | C ₃ H ₆ O ₃ | Conc. 90% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 | |
| Nitric acid | HNO ₃ | Conc. 65% | 2min | A5 | B1 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A6 | B2 | C1 | D0 | E0 | F5 | G0 | |
| | | | 24 h | A6 | B2 | C1 | D1 | E1 | F6 | G0 | |
| Nitric acid | HNO ₃ | 30% | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A5 | B1 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A5 | B2 | C1 | D0 | E0 | F0 | G0 | |
| Oxalic acid | C ₂ H ₂ O ₄ | 10% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Sulphuric acid | H ₂ SO ₄ | Conc. 98% | 2min | A5 | B1 | C0 | D0 | E0 | F0 | G0 | H* |
| | | | 1 h | A6 | B2 | C1 | D1 | E1 | F5 | G0 | H* |
| | | | 24 h | A6 | B2 | C1 | D2 | E1 | F5 | G0 | H* |
| Sulphuric acid | H ₂ SO ₄ | 30% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 | |
| Organic solvents | | | | | | | | | | | |
| Acetone | C ₃ H ₆ O | | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B1 | C1 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B1 | C1 | D1 | E5 | F5 | G0 | |
| Acetonitrile | CH ₃ CN | | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E5 | F5 | G0 | |
| Carbon tetrachloride | CCl ₄ | | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B0 | C0 | D0 | E0 | F1* | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F1* | G0 | |
| Chloroform | CHCl ₃ | | 2min | A0 | B1 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B2 | C1 | D0 | E5 | F0 | G0 | |
| | | | 24 h | A0 | B2 | C1 | D1 | E5 | F5 | G0 | |
| Cyclohexane | C ₆ H ₁₂ | | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Cyclohexanone | C ₆ H ₁₀ O | | 2 min | A0 | B1 | C1 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B2 | C2 | D3 | E5 | F0 | G0 | H |
| | | | 24 h | A0 | B2 | C2 | D3 | E6 | F6 | G0 | H |
| Dichloroethylene | C ₂ H ₂ Cl ₂ | | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E5 | F5 | G0 | |
| Methylene Chloride | CH ₂ Cl ₂ | | 2min | A0 | B1 | C1 | D0 | E0 | F0 | G0 | H* |
| | | | 1 h | A0 | B2 | C1 | D0 | E5 | F5 | G0 | H* |
| | | | 24 h | A0 | B2 | C1 | D0 | E5 | F5 | G0 | H* |
| Ethanol | C ₂ H ₅ OH | | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 | |

| Organic solvents,cont. | | | | | | | | | | | |
|---------------------------------|---|------------------------|-------|----|----|----|----|----|-----|----|----|
| Ethyl acetate | C ₄ H ₈ O ₂ | | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B1 | C1 | D0 | E5 | F1* | G0 | |
| | | | 24 h | A0 | B1 | C1 | D0 | E5 | F5 | G0 | |
| Ethylene glycol | C ₂ H ₆ O ₂ | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Diethyl ether | (C ₂ H ₅) ₂ O | | 2 min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B0 | C0 | D0 | E0 | F1* | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F5 | G0 | |
| n-Hexane | C ₆ H ₁₄ | | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Formaldehydesolution | CH ₂ O | 37 % | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Methanol | CH ₃ OH | | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 | |
| Methyl ethyl ketone | C ₄ H ₈ O | | 2 min | A0 | B1 | C1 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B2 | C2 | D2 | E5 | F0 | G0 | H* |
| | | | 24 h | A0 | B2 | C2 | D3 | E5 | F5 | G0 | H* |
| Pet.ether (Ligroin) 80-110°C | CAS-nr: 8032-32-4 | | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Tetrachloroethylene | C ₂ Cl ₄ | | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B0 | C0 | D0 | E0 | F1* | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E5 | F1* | G0 | |
| Toluene | C ₇ H ₈ | | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B1 | C0 | D0 | E0 | F1* | G0 | |
| | | | 24 h | A0 | B2 | C1 | D0 | E5 | F5 | G0 | |
| Trichlorethylene | C ₂ HCl ₃ | | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B1 | C1 | D0 | E5 | F1* | G0 | |
| | | | 24 h | A0 | B1 | C1 | D0 | E5 | F5 | G0 | |
| White spirit | EG/EC/EF-no: 265-191-7 | | 2min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Xylene | C ₈ H ₁₀ | | 2 min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B1 | C1 | D0 | E5 | F5 | G0 | |
| Alkali (Bases) | | | | | | | | | | | |
| Ammonia solution | NH ₃ | 25% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Calcium hydroxide | Ca(OH) ₂ | 10% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Sodium hydroxide | NaOH | 50% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Sodium hydroxide | NaOH | 10% | 1 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B2 | C0 | D0 | E0 | F0 | G0 | H |
| Salt solutions | | | | | | | | | | | |
| Ammonium carbonate | (NH ₄) ₂ CO ₃ | 10% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Ammonium iron (III) sulphate | NH ₄ Fe(SO ₄) ₂ | 10% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Calcium Chloride | CaCl ₂ | Saturated | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Cobaltous chloride | CoCl ₂ | 10% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Copper (II) sulphate | CuSO ₄ | 10% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Ferrous (II) chloride | FeCl ₂ | 10% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A1 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Ferric (III) chloride | FeCl ₃ | 10% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A1 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Potassium iodide | KI | 10% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Potassium oxalate | K ₂ C ₂ O ₄ | 10% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Potassium permanagnate | KMnO ₄ | 5% in H ₂ O | 2 min | A5 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A6 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Silver nitrate | AgNO ₃ | 2% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A5 | B0 | C0 | D0 | E0 | F0 | G0 | |

| Salt solutions cont. | | | | | | | | | | |
|---|---|-----------------|------|----|----|----|----|----|----|----|
| Sodium carbonate | Na ₂ CO ₃ | 20% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Sodium thiosulphate | Na ₂ S ₂ O ₃ | 10% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 |
| Sodium sulphite | Na ₂ SO ₃ | 10% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Medical chemicals | | | | | | | | | | |
| Aniline blue | | 2,5% in ethanol | 1 h | A6 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | | 24 h | A6 | B0 | C0 | D0 | E0 | F0 | G0 |
| Betadine skin cleanser | | 75mg/ml | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | | 24 h | A5 | B0 | C0 | D0 | E0 | F0 | G0 |
| Bromocresol green | | 0,04 % | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Eosin | | 1 % in Ethanol | 1 h | A6 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | | 24 h | A6 | B0 | C0 | D0 | E0 | F0 | G0 |
| Glutaraldehyde | | 25% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Hematoxylin | | 5% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | | 24 h | A5 | B0 | C0 | D0 | E0 | F0 | G0 |
| Hibitane | | 0,5% | 1 h | A5 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | | 24 h | A6 | B0 | C0 | D0 | E0 | F0 | G0 |
| Iodine | I ₂ | 2 % in ethanol | 2min | A6 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | | 1 h | A6 | B0 | C0 | D0 | E0 | F0 | G0 |
| Iodoform | | 1% in ethanol | 1 h | A6 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | | 24 h | A6 | B0 | C0 | D0 | E0 | F0 | G0 |
| Methylrosanilinium | | 0,1% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | | 24 h | A6 | B0 | C0 | D0 | E0 | F0 | G0 |
| Disinfectants/cleaning compounds | | | | | | | | | | |
| Product | Manuf./Rep. | | | | | | | | | |
| Buraton 10F | Schülke & Mayr | 1% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 10% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Citrosteril | Fresenius | Conc. | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Debisan | Nordex | 1 % | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 10% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Decon-Spore 200 Plus | Veltek Associates, Inc | 0,5 % | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | 5 % | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Dialox | Gambro | Conc. | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 |
| Gevisol | Schülke & Mayr | 0,5% | 24 h | A5 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 5% | 24 h | A6 | B0 | C0 | D0 | E0 | F0 | G0 |
| Incidur | Henkel | 0,5% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 3% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Lycetol AF | Schülke & Mayr | 1% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 5% | 24 h | A5 | B0 | C0 | D0 | E0 | F0 | G0 |
| Melsept | B Braun | 1% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 5% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Perform | Schülke & Mayr | 0,75% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 2,5% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Sekumatic | Henkel | 0,5% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 5% | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 |
| Sekusept Plus | Henkel | 1% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 5% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Spitacid | Henkel | Conc. | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 |
| Terralin N | Schülke & Mayr | 1% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 10% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 10% | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 |
| Tiutol KF | B. Braun | 3% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 10% | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 |
| Virkon S | Sterisol AB | 1% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| ''- | | 2,5% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| Incidin Plus | Ecolab | 1% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |
| | | 5% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 |

| Disinfectants/cleaning compounds. cont. | | | | | | | | | | | |
|---|---|------|-------|----|----|----|----|----|----|----|--|
| Product | Manuf./Rep. | | | | | | | | | | |
| Incidin Extra N | Ecolab | 1% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | 5% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Mikrobac forte | BODE Chemi | 1% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | 5% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Hexaquart plus | B. Braun | 1% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | 2,5% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Miscellaneous chemicals | | | | | | | | | | | |
| EDTA | C ₁₀ H ₁₆ N ₂ O ₈ | 10% | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Glycerol | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Hydrogen peroxide | H ₂ O ₂ | 30% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Olive oil | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| Phenol | C ₆ H ₆ O | 5% | 2 min | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 1 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 | |
| Sodium hypochlorite | NaOCl | 12% | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B1 | C0 | D0 | E0 | F0 | G0 | |
| Brake fluid Super DOT 4 | APE Components AB | Conc | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D1 | E5 | F0 | G0 | |
| Hydraulic fluid DET 26 | Mobil | Conc | 1 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |
| | | | 24 h | A0 | B0 | C0 | D0 | E0 | F0 | G0 | |

*The swelling disappears after 1-2 days.

H* Slight damage to polyurethane surface.

H Total damage to polyurethane surface.

RESISTANCE TO CHEMICALS

Key

- A0 No change in lightness or colour
- A1 Somewhat lighter surface
- A2 Lighter surface
- A3 Somewhat darker surface
- A4 Darker surface
- A5 Somewhat discoloured surface
- A6 Discoloured surface

- B0 No change in gloss or matness
- B1 Somewhat mat surface
- B2 Mat surface
- B3 Somewhat glossy surface
- B4 Glossy surface

- C0 No change in patchiness
- C1 Somewhat patchy or spotty surface
- C2 Patchy or spotty surface

- D0 No change in evenness
- D1 Somewhat uneven or porous surface
- D2 Uneven and porous surface
- D3 Somewhat crackled surface
- D4 Crackled surface

- E0 No brittleness, stickiness or softening
- E1 Some surface brittleness
- E2 Brittle surface
- E3 Some surface stickiness
- E4 Sticky surface
- E5 Somewhat softened
- E6 Softened

- F0 No change in size or flatness
- F1 Some swelling
- F2 Swelling
- F3 Slight shrinkage
- F4 Shrinkage
- F5 Some bulging
- F6 Bulging or distortion

- G0 No delamination
- G1 Delamination of two or more layers

- H Other changes noted (text en clair)